

Capture and Review Patient Photos and Media




Audience: physicians, mid-level providers, nurses, and other clinical support staff


Below are the steps to capture images in Haiku/Canto

Capture clinical images


1. Select a patient from the **Schedule/Patient List**.



2. Tap on the patient's picture at the top of the screen.



3. In the Media Capture activity, tap **Capture Clinical Media** and use the camera on your device to take a picture.



4. After taking a picture. Tap **Use Photo**.




Capture and Review Patient Photos and Media




5. Select a Document Type, if necessary, add a Description. When complete you can add another picture or save.

6. Tap **Save**. You can review media captures in the Media Manager or the Chart Review activities in Hyperspace.



Review clinical images and media

1. In the patient's chart, tap .
2. Scroll through the list of media files or, to find a specific media file, search using information such as:
 - Document type
 - The user who imported the file
 - Keywords from the image's description
 - The date the file was imported
3. Tap a photo or document to see a full-screen view, where you can zoom in and out or pan around the image.