

Dragon: Sound Settings for PowerMic and Video Visits




Audience: All Providers


Set your computer sound for Dragon PowerMic and video visits.

Sound Settings for PowerMic and Video Visits

1. Hit Windows **Start** button. Go to **Settings**.




2. Click **System**.



Dragon: Sound Settings for PowerMic and Video Visits


3. Go to Sound.

- **Sound:**
 - Set to **Speaker/HP (Realtek High Definition)** to hear your patient's voice from your computer audio.
- **Input:**
 - Set to **Microphone – PowerMic-NS** so you can:
 - Dictate in Dragon.
 - Speak to your patient during the video visit.



Dragon: Sound Settings for PowerMic and Video Visits

4. Make sure you **Log Off** from Dragon:
 - Click the **flame**, then click **Log Off**.



5. The **splash screen** opens.
6. Choose your **Microphone** source as **PowerMic-NS** (handheld) or **Nuance PowerMic Mobile** (cell phone). **Note:** Make sure computer audio (**Realtek**) is **not** selected.

