

Amnioinfusion [3386]

Nursing

Vitals

Vital signs - T/P/R/BP Routine, Per unit protocol, L&D Pre-Delivery

Activity

- Bed rest Routine, Until discontinued, Starting S
Bathroom Privileges:
L&D Pre-Delivery
- Bed rest with bathroom privileges Routine, Until discontinued, Starting S
Bathroom Privileges: with bathroom privileges
L&D Pre-Delivery
- Bed rest with bathroom privileges Routine, Until discontinued, Starting S
Bathroom Privileges: with bathroom privileges
For Bowel Movement Only, L&D Pre-Delivery

Nursing

- Begin Amnioinfusion Routine, Until discontinued, Starting S, L&D Pre-Delivery
- Monitor fetal heart tones Routine, Once
Type: Continuous
L&D Pre-Delivery
- IUPC Monitor Routine, Continuous, L&D Pre-Delivery
- Strict intake and output Routine, Every 2 hours
Strict for Amniotic fluid, L&D Pre-Delivery
- Reposition patient to aid in output Routine, Now then every 2 hours, L&D Pre-Delivery
- Weigh Towel/Chuck Add-On, As needed
Every Peri-Pad Change, L&D Pre-Delivery
- Discontinue Amnioinfusion Routine, Until discontinued, Starting S
Discontinue Amnioinfusion if no output noted, L&D
Pre-Delivery

Notify

- Notify Admitting MD No amniotic fluid output Routine, Until discontinued, Starting S, L&D Pre-Delivery
- Notify Admitting MD fluid retention Routine, Until discontinued, Starting S, L&D Pre-Delivery

Intrauterine Fluids

Intrauterine Fluids (Single Response)

- () Lactated Ringers Bolus and Infusion **"Followed by" Linked Panel**
- lactated Ringer's intrauterine bolus from bag 500 mL, intrauterine, for 30 Minutes, once, For 1 Doses, L&D Pre-Delivery
Do not administer through the IV Pump
- lactated ringer's intrauterine infusion 125 mL/hr, intrauterine, continuous, Starting H+30 Minutes, L&D
Pre-Delivery
Do not administer through the IV Pump
- () Sodium Chloride 0.9% Bolus and Infusion **"Followed by" Linked Panel**
- sodium chloride 0.9 % intrauterine bolus from bag 500 mL, intrauterine, for 30 Minutes, once, For 1 Doses, L&D Pre-Delivery
Do not administer through the IV Pump
- sodium chloride 0.9 % intrauterine infusion 125 mL/hr, intrauterine, continuous, Starting H+30 Minutes, L&D
Pre-Delivery
Do not administer through the IV Pump